

Perform Sterile Alcohol EP

STERIL, ALKOHOLBASERET OPLØSNING

- Hurtig, effektiv og afsætter ingen rester
- Steril (gamma steriliseret)
microbial filtreret (0,2 µm)
- Renrum grade A/B, ISO 5/6
- Klar til brug i handy og steril 500 ml sprayflaske
- Farve og parfumefri
- Effektivitetstest efter europæisk standard

MIKROBIOLOGISK EFFEKTIVITET

- Bactericidal EN1276, EN13697
Kontaktid: 1 min.
- Fungicidal EN1650, EN13697
Kontaktid: 2 min.
- Verucidal
Kontaktid: 30 sek.



VARENUMMER	TYPE	PAKNING
121802	500 ml sprayflaske	10 x 500 ml/karton



Sterile, alcohol-based ready to use solution for disinfecting surfaces in cleanrooms.

perform® sterile alcohol EP

Our Plus

- Fast, effective and residue free
- Sterile (γ -irradiated) and microbial filtrated (0.2 μm)
- Suitable for cleanrooms grade A/B, ISO class 5/6
- Ready to use in a handy, sterile 500 ml spray system
- Colour and perfume free
- Efficacy tested according to European norms (EN)

Application areas

Disinfection of surfaces and equipment/devices in cleanrooms. Especially for Grade A and B clean areas, according to:

- EU Guide to Good Manufacturing Practice for medicinal products Annex 1: Manufacture of sterile medicinal products, March 2009.
- FDA Draft Guidance for Industry: Sterile Drug Products produced by aseptic processing- current Good Manufacturing practice, U.S. Department of Health and Human Services; Food and Drug Administration, September 2004.

Instructions for use

Apply perform® sterile alcohol EP undiluted to objects and surfaces and allow to take effect. Ensure complete wetting. The product dries off quickly. perform® sterile alcohol EP can be used for glove disinfection. Maximum application rate 50 ml/sq.m. Do not use to treat surfaces that are particularly sensitive to alcohol, such as acrylic glass.

Efficacy	Concentration	Contact time
bactericidal EN1276, EN13697	ready-to-use	1 min.
fungicidal EN1650, EN13697	ready-to-use	2 min.
virucidal against enveloped viruses in accordance with DVG (German Registered Association for Combating Viral Diseases)//RKI Guideline	ready-to-use	30 sec.



perform® sterile alcohol EP

Product data

Composition:

100 g solution contains the following active substances: 25 g Ethanol (94 %), 35 g Propan-1-ol.

Chemical-physical data

Color	colourless
Density	ca. 0,89 g/cm ³ / 20 °C
Flash point	27 °C / Method : DIN 51755 Part 1
Form	liquid
pH	ca. 7 / 20 °C

Special advice

Use disinfectants safely. Always read the label and product information before use.

The sterile 500 ml spray system offers the following advantages:

- convenient size for storage
- ergonomically shaped spray pump for ease of use
- remains sterile after opening

The bottles are sterile (with irradiation indicator) and doublewrapped individually. perform® sterile alcohol EP comprises alcohol of pharmaceutical quality in accordance with the Ph. Eur.

Information for order

Item	Delivery form	Item no.
Perform® sterile alcohol EP 500 ml bottle	10/Carton	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Application aids

Application aids	Item no.
Perform® spray bottle holder	693540
Perform® spray bottle holder- plate	693541



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.



A company of the Air Liquide-Group

Related Products

- Perform® sterile dry wipes
- Perform® sterile mix dry wipes
- Perform® sterile premium dry wipes

Environmental information

Schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: www.schuelke.com.

For individual questions:

Customer Sales Service

Phone: +49 40 52100-666

E-Mail: info@schuelke.com

Schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40- 52100- 0
Fax +49 40- 52100- 318
www.schuelke.com
info@schuelke.com